



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz separatory funnel, cylindrical, with solid glass plug, 100 ml, NS 29/32, bore 2.5 mm, graduation 5 ml

**€103.28
plus VAT &
Shipping**

Product Images



Description

Cylindrical shape with smooth drain tube, with NS sleeve and plastic stopper (PE).

With solid NS glass plug, with threaded stopcock safety device, with graduation for reading off the Zirka volume.

Material: Borosilicate 3.3

Additional Information

| | |
|------------------------|--|
| No. | LZ-0839 |
| Manufacturer (Brand) | Lenz |
| old neoLab article no. | B-2839 |
| EAN | 4058072305796 |
| Transport temperature | Room temperature |
| autoclavable at 121°C | Yes |
| Color | transparent (clear) |
| Material | Borosilicate glass |
| Volume | 100 ml |
| sterile | No |
| Form | cylindrical |
| Type closure/ lid | NS plug |
| Funnel features | with closure/lid with standard ground joint with NS stopcock/valve graduated |
| Color graduation | white |
| Material chucks | Borosilicate glass |
| Standard ground | 29/32 |
| Division | 5 ml |
| Material closure/ lid | Polyethylene (PE) |

