




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz separatory funnel, cylindrical, with PTFE
plug, 50ml, NS 29/32, bore 2.5 mm, graduation 2
ml**

**€130.22
plus VAT &
Shipping**



Description

Cylindrical shape with smooth drain pipe, with NS sleeves and plastic plugs (PE), with PTFE plug.
With threaded stopcock safety device, with graduation for reading off the Zirka volume.
Material: borosilicate glass 3.3

Additional Information

No.	LZ-0843
Manufacturer (Brand)	Lenz
old neoLab article no.	B-2843
EAN	4058072305833
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Volume	50 ml
sterile	No
Form	cylindrical
Type closure/ lid	NS plug
Funnel features	with closure/lid with standard ground joint with NS stopcock/valve graduated
Color graduation	white
Material chucks	Polytetrafluoroethylene (PTFE)
Standard ground	29/32
Division	2 ml
Material closure/ lid	Polyethylene (PE)

