



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**Lenz separatory funnels, cylindrical, with PTFE
plug, 100 ml, NS 29/32, bore 2.5 mm, graduation
5 ml**

**€123.16
plus VAT &
Shipping**



Description

Cylindrical shape with smooth drain pipe, with NS sleeves and plastic plugs (PE), with PTFE plug.
With threaded stopcock safety device, with graduation for reading off the Zirka volume.
Material: borosilicate glass 3.3

Additional Information

No.	LZ-0844
Manufacturer (Brand)	Lenz
old neoLab article no.	B-2844
EAN	4058072305840
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Volume	100 ml
sterile	No
Form	cylindrical
Type closure/ lid	NS plug
Funnel features	with closure/lid with standard ground joint with NS stopcock/valve graduated
Color graduation	white
Material chucks	Polytetrafluoroethylene (PTFE)
Standard ground	29/32
Division	5 ml
Material closure/ lid	Polyethylene (PE)

