



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**Lenz separatory funnels, cylindrical, with PTFE  
plug, 1,000 ml, NS 29/32, bore 6.0 mm,  
graduation 20 ml**

**€239.18  
plus VAT &  
Shipping**

## Product Images

---



## Description

---

Cylindrical shape with smooth drain pipe, with NS sleeves and plastic plugs (PE), with PTFE plug.  
With threaded stopcock safety device, with graduation for reading off the Zirka volume.  
Material: borosilicate glass 3.3

## Additional Information

---

No.	LZ-0847
Manufacturer (Brand)	Lenz
old neoLab article no.	B-2847
EAN	4058072305871
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Volume	1000 ml
sterile	No
Form	cylindrical
Type closure/ lid	NS plug
Funnel features	with closure/lid with standard ground joint with NS stopcock/valve graduated
Color graduation	white
Material chucks	Polytetrafluoroethylene (PTFE)
Standard ground	29/32
Division	20 ml
Material closure/ lid	Polyethylene (PE)

