



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



---

## Lenz dropping funnel, conical, with solid glass plug, 100 ml, NS 19/26, bore 2.5 mm

**€66.42**  
**plus VAT &**  
**Shipping**

## Product Images

---



## Description

---

Dropping funnel, conical, with glass plug

Conical shape with NS sleeve and NS core, with solid NS glass plug.

With threaded stopcock safety device, with labeling field, without graduation.

Material: borosilicate glass 3.3.

## Additional Information

---

No.	LZ-0893
Manufacturer (Brand)	Lenz
old neoLab article no.	B-2893
EAN	4058072306335
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Volume	100 ml
Form	conical
Funnel features	with standard ground joint
Standard ground sleeve	19/26
Standard ground core	19/26

