



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz dropping funnel cylindrical, with solid glass plug, 50 ml, NS 29/32, bore 2.5 mm, graduation 1 ml

**€134.02
plus VAT &
Shipping**

Product Images



Description

Dropping funnel, cylindrical, with glass plug.

Cylindrical shape with NS sleeves and NS core, with solid NS glass plug, with threaded stopcock safety device.

With graduation for reading off the Zirka volume, with labeling field.

Material: borosilicate glass 3.3.

Additional Information

No.	LZ-0862
Manufacturer (Brand)	Lenz
old neoLab article no.	B-2862
EAN	4058072306021
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Volume	50 ml
Form	cylindrical
Funnel features	with standard ground joint
Standard ground sleeve	29/32
Standard ground core	29/32

