



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e


Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**Lenz dropping funnel cylindrical, with solid glass plug, 1.000 ml, NS 29/32, bore 6,0 mm, graduation 20 ml**

**€254.10  
plus VAT &  
Shipping**

## Product Images



## Description

---

Dropping funnel, cylindrical, with glass plug.

Cylindrical shape with NS sleeves and NS core, with solid NS glass plug, with threaded stopcock safety device. With graduation for reading off the Zirka volume, with labeling field.

Material: borosilicate glass 3.3.

## Additional Information

---

No.	LZ-0867
Manufacturer (Brand)	Lenz
old neoLab article no.	B-2867
EAN	4058072306076
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Volume	1000 ml
Form	cylindrical
Funnel features	with standard ground joint
Standard ground sleeve	29/32
Standard ground core	29/32

