



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Lenz dropping funnel cylindrical, with pressure equalization tube, PTFE plug, 250 ml, NS 29/32, bore 4.0 mm, graduation 5 ml

**€223.20
plus VAT &
Shipping**

Product Images



Description

Dropping funnel, cylindrical, with PTFE plug, with pressure compensation tube

Cylindrical shape with NS sleeve and NS core, with pressure compensation tube.

With PTFE plug, with threaded stopcock safety device, with graduation for reading off the Zirka volume, with inscription field.

Material: borosilicate glass 3.3.

Additional Information

No.	LZ-0889
Manufacturer (Brand)	Lenz
old neoLab article no.	B-2889
EAN	4058072306298
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Borosilicate glass
Volume	250 ml
Form	cylindrical
Funnel features	with standard ground joint
Standard ground sleeve	29/32
Standard ground core	29/32

