



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® Storage bottles for electrodes, 22 x 68 mm high, 5 pcs/pack

**€16.50**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Many pH electrodes, must not dry out and are usually stored in 3 molar KCL solution. These storage bottles for electrodes hold approx. 20 ml of solution. The silicone seal holds the electrode firmly and stably in the bottle. Simply clamp the bottle into a clamp and use it to store the electrode without worry.  
5 bottles/pack; without electrolyte solution.

## Additional Information

---

No.	3-1951
Manufacturer (Brand)	neolab
EAN	4058072121143
Transport temperature	Room temperature
Color	transparent (clear)
Volume electrolyte solution	20 ml

