



neoLab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neoCulture cell culture flasks 25 cm², 65 ml,
sterile, PS, vent cap, box of 200 (20x10pcs.)**

Special Price

€118.00 was

€153.70

**plus VAT &
Shipping**

Product Images



Description

The cell and tissue culture flasks with treated surface are perfect for adherent cell growth.

The surface is flat and free of streaking to maximize the usable growth area.

The standard closed caps can be used in an open system.

An innovative angled neck design provides good access when working with pipettes and cell scrapers. The bottles can be easily stacked, and a special outer area near the neck is used for marking. Graduation is engraved on both sides of the bottle.

Sterilized by gamma radiation

DNase-/RNase-free

pyrogen-free

Additional Information

No.	C-8200
Manufacturer (Brand)	neoCulture
EAN	4058072154660
Transport temperature	Room temperature
Material	Polystyrene (PS)
Volume	50 ml
sterile	Yes
Type closure/ lid	Screw cap with ventilation
Color closure/ lid	white
Fluid behavior	hydrophobic
Material closure/ lid	Polyethylene (PE)
Closure / lid properties	closed
Growth area	25.00 cm ²

