



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neoCulture cell culture flasks 75 cm², 293 ml, PS, vent cap, box of 100 (20x5) sterile

Special Price

€128.00 was

€166.10

**plus VAT &
Shipping**

Product Images



Description

The cell and tissue culture flasks with treated surface are perfect for adherent cell growth.
The surface is flat and free of streaking to maximize the usable growth area.
The standard closed caps can be used in an open system.
An innovative angled neck design provides good access when working with pipettes and cell scrapers. The bottles can be easily stacked, and a special outer area near the neck is used for marking. Graduation is engraved on both sides of the bottle.
Sterilized by gamma radiation
DNase-/RNase-free
pyrogen-free

Additional Information

No.	C-8201
Manufacturer (Brand)	neoCulture
EAN	4058072155216
Transport temperature	Room temperature
Material	Polystyrene (PS)
Volume	250 ml
sterile	Yes
Type closure/ lid	Screw cap with ventilation
Fluid behavior	hydrophobic
Closure / lid properties	closed
Growth area	75.00 cm2
Color closure/ lid	white
Material closure/ lid	Polyethylene (PE)

