



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neoCulture cell culture flasks 175 cm², 580 ml,
PS, gas exchange cap, box of 40 (8x5) sterile**

Special Price
€84.70 was
€110.00
plus VAT &
Shipping

Product Images



Description

The cell and tissue culture flasks with treated surface are perfect for adherent cell growth.

The surface is flat and free of streaking to maximize the usable growth area.

The vented caps can be used in a closed system.

An innovative angled neck design provides good access when working with pipettes and cell scrapers.

The bottles can be easily stacked, and a special outer area near the neck is used for marking. Graduation is engraved on both sides of the bottle.

Sterilized by gamma radiation

DNase/RNase free

Non-pyrogenic

Additional Information

No.	C-8205
Manufacturer (Brand)	neoCulture
EAN	4058072154790
Transport temperature	Room temperature
Material	Polystyrene (PS)
Volume	580 ml
sterile	Yes
Type closure/ lid	Screw cap with gas exchange
Fluid behavior	hydrophobic
Closure / lid properties	closed
Growth area	175.00 cm ²
Color closure/ lid	white
Material closure/ lid	Polyethylene (PE)

