



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Narrow neck screw cap bottle 250 ml, clear glass, DIN 22, without screw cap

€2.09 As low as

€1.75

**plus VAT &
Shipping**



Description

Narrow-neck screw-top vials made of clear soda glass, round, according to DIN.
Fits perfectly into the neoLab shipping boxes.

Accessories: Screw thread caps according to DIN (Order No. 2-1490 et seq.)

Additional Information

No.	DT-0068
Manufacturer (Brand)	neoLab
old neoLab article no.	2-1453
EAN	4058072228453
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Soda-lime glass (normal glass) 2
Volume	250 ml
Color bottle	transparent (clear)
Height	148.5 mm
Material	Soda-lime glass (normal glass) 2
Thread	DIN 22
Neck size	narrow
DM bottle	62 mm
DM Neck	22 mm
Material bottle	Soda-lime glass (normal glass) 2

