



neolab Migge GmbH

Rischerstr. 7-9

69123 Heidelberg

Deutschland

+49 (0)6221 / 8442-44

<https://www.neolab.de>

Umsatzsteuer-Identifikationsnummer:

DE 143 450 657



neolab® Narrow neck screw cap bottle 500 ml, clear glass, DIN 25, without screw cap

€3.17 As low as

€2.66

plus VAT & Shipping



Description

Narrow-neck screw-top vials made of clear soda glass, round, according to DIN.

Fits perfectly into the neoLab shipping boxes.

Accessories: Screw thread caps according to DIN (Order No. 2-1490 et seq.)

Additional Information

No.	DT-0069
Manufacturer (Brand)	neoLab
old neoLab article no.	2-1454
EAN	4058072228460
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Soda-lime glass (normal glass) 2
Volume	500 ml
Color bottle	transparent (clear)
Height	183.7 mm
Material	Soda-lime glass (normal glass) 2
Thread	DIN 25
Neck size	narrow
DM bottle	76 mm
DM Neck	25 mm
Material bottle	Soda-lime glass (normal glass) 2

