



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® Narrow neck screw cap bottle 1000 ml,
clear glass, DIN 28, without screw cap**

€4.68 As low as
€3.77
**plus VAT &
Shipping**

Product Images



Description

Narrow-neck screw-top vials made of clear soda glass, round, according to DIN.
Fits perfectly into the neoLab shipping boxes.

Accessories: Screw thread caps according to DIN (Order No. 2-1490 et seq.)

Additional Information

No.	DT-0070
Manufacturer (Brand)	neoLab
old neoLab article no.	2-1455
EAN	4058072228477
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Material	Soda-lime glass (normal glass) 2
Volume	1000 ml
Color bottle	transparent (clear)
Height	218.3 mm
Material	Soda-lime glass (normal glass) 2
Thread	DIN 28
Neck size	narrow
DM bottle	97 mm
DM Neck	28 mm
Material bottle	Soda-lime glass (normal glass) 2

