



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® PET narrow neck bottles glass cl. 200ml

€30.00
plus VAT &
Shipping



Description

Plastic bottles made of crystal-clear polyester.
Polypropyleneclosures.

Food safe, high chemical resistance.
Largely diffusion resistant.
Temperature range: -40°C to +60°C

Additional Information

No.	3-2080
Manufacturer (Brand)	neoLab
EAN	4058072018832
Transport temperature	Room temperature
Color	transparent (crystal clear)
Material	Polyethylene terephthalate (PET)
Volume	200 ml
sterile	No
Color bottle	transparent (crystal clear)
Height	120 mm
Type closure/ lid	Screw cap
TBST MAX	60 °C
TBST MIN	-40 °C
Bottles properties	with closure/lid closure autoclavable
Material	Polyethylene terephthalate (PET)
Thread	DIN 28
Neck size	narrow
Closure / lid properties	closed
DM bottle	55 mm
DM Neck	22 mm
Color closure/ lid	white
Material bottle	Polyethylene terephthalate (PET)
Material closure/ lid	Polypropylene (PP)

