



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® Narrow neck screw cap bottle brown, 100 ml without screw cap DIN 22

**€1.15** As low as  
**€1.08**  
plus VAT &  
Shipping

## Product Images

---



## Description

---

Narrow neck screw top vials made of brown soda glass, round, acc. to DIN.

Fits perfectly into the neoLab shipping boxes.

Accessories: Screw thread caps acc. to DIN (Order No. 2-1490 et seq.)

## Additional Information

---

No.	DT-0071
Manufacturer (Brand)	neoLab
old neoLab article no.	2-1462
EAN	4058072228484
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	brown
Material	Soda-lime glass (normal glass) 2
Volume	100 ml
Color bottle	brown
Height	117 mm
Material	Soda-lime glass (normal glass) 2
Thread	DIN 22
Neck size	narrow
DM bottle	45 mm
DM Neck	22 mm
Material bottle	Soda-lime glass (normal glass) 2

