



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Narrow neck screw top bottle brown, 1000 ml without screw cap DIN 28

€4.43 As low as

€3.71

**plus VAT &
Shipping**

Product Images



Description

Narrow neck screw top vials made of brown soda glass, round, acc. to DIN.
Fits perfectly into the neoLab shipping boxes.
Accessories: Screw thread caps acc. to DIN (Order No. 2-1490 et seq.)

Additional Information

No.	DT-0074
Manufacturer (Brand)	neoLab
old neoLab article no.	2-1465
EAN	4058072228514
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	brown
Material	Soda-lime glass (normal glass) 2
Volume	1000 ml
Color bottle	brown
Height	218.3 mm
Material	Soda-lime glass (normal glass) 2
DM bottle	97 mm
DM Neck	28 mm
Thread	DIN 28
Neck size	narrow
Material bottle	Soda-lime glass (normal glass) 2

