



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



DURAN® pressure plus laboratory bottle with DIN thread GL 45, 1000 ml vacuum/pressure- tested up to 1.5 bar

€26.81 As low as
€20.62
plus VAT &
Shipping

Product Images



Description

Duran® pressure plus laboratory bottle with DIN thread GL 45.

With easy-to-read graduation and large text field. Burnt-in, therefore very long-lasting pressure image.

Tested pressure resistance according to DIN ISO 1595 (GS) . Pressure resistance from -1 to +1.5 bar (at max. 140°C) due to modified geometry. With pressure load thermal shock resistance 30 K and max. +140°C service temperature.

Blue graduation for visual differentiation from standard bottles.

Capacity: 1000 ml

Diameter: 101 mm

Height : 225 mm

Additional Information

No.	SD-0892
Manufacturer (Brand)	DURAN
old neoLab article no.	2-3050
EAN	4058072346119
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Volume	1000 ml
Color bottle	transparent (clear)
Height	225 mm
TBST MAX	140 °C
TBST MIN	-140 °C
Bottles properties	with high hydrolytic resistance with special thermal shock resistance with special pressure resistance with batch certificate with labeling field graduated
Thread	GL 45
Neck size	narrow
DM bottle	101 mm
DM Neck	45 mm
Material bottle	Borosilicate glass (1st hydrolytic class)

