



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Vitlab laboratory bottle, PP with tamper evident closure, PP, GL 45, 500 ml, 6 pcs/pack

**€45.00**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

The all-round bottle made of polypropylene for the storage of liquids in the laboratory. Due to their innovative design, the bottles have an excellent grip and thus make work easier compared to conventional laboratory bottles. Due to optimized sealing system of bottle thread and screw cap leakage safety is guaranteed. Including GL45 tamper-evident cap **Volume: 500 ml**

## Additional Information

No.	LV-0092
Manufacturer (Brand)	Vitlab
old neoLab article no.	2-0124
EAN	4058072289157
Transport temperature	Room temperature
Color	transparent (milky)
Material	Polypropylene (PP)
Volume	500 ml
sterile	No
Color bottle	transparent (milky)
Height	192 mm
Type closure/ lid	Tamper-evident closure
TBST MAX	125 °C
TBST MIN	0 °C
Bottles properties	with closure/lid graduated closure autoclavable bottle autoclavable
Material	Polypropylene (PP)
Thread	GL 45
Neck size	narrow
Closure / lid properties	closed
DM bottle	77 mm
DM Neck	45 mm
Spare part	LV-0975
Color closure/ lid	red
Material bottle	Polypropylene (PP)
Material closure/ lid	Polypropylene (PP)

