




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Laboratory bottles amber glass borosilicate 1000 ml, 10 pcs./pack

€235.60
plus VAT &
Shipping

Product Images



Description

SIMAX® glass with brown coating,DIN thread, with screw cap and pouring ring made of PP/blue (max. 140°C) Batch identification (Retrace Code) for traceability autoclavable!

Additional Information

No.	E-1824
Manufacturer (Brand)	neoLab
EAN	4058072067335
Transport temperature	Room temperature
Color	brown
Material	Borosilicate glass
Volume	1000 ml
sterile	No
Color bottle	brown
Height	225 mm
Type closure/ lid	Screw cap
Bottles properties	with closure/lid with batch certificate with pouring ring graduated closure autoclavable bottle autoclavable
Material	Borosilicate glass
DM bottle	101 mm
DM Neck	45 mm
Color closure/ lid	blue
Thread	GL 45
Neck size	narrow
Material bottle	Borosilicate glass
Material closure/ lid	Polypropylene (PP)
Closure / lid properties	closed

