



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



neolab® Laboratory bottles oh. Cap 250 ml ISO 4796 Boro glass 3.3 GL 45 PU 10 pieces

€23.90
plus VAT &
Shipping

Product Images



Description

Graduated laboratory bottles with screw thread made of borosilicate glass 3.3 (acc. to ISO 4796 DIN 168 part 1).

Additional Information

No.	2-3061
Manufacturer (Brand)	neoLab
EAN	4058072000479
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Volume	250 ml
Color bottle	transparent (clear)
Height	143 mm
Bottles properties	with high hydrolytic resistance with labeling field with pouring ring graduated
Thread	DIN 45
Neck size	narrow
DM bottle	70 mm
DM Neck	45 mm
Material bottle	Borosilicate glass (1st hydrolytic class)

