



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## neolab® Laboratory bottles oh. Cap 500 ml ISO 4796 Boro glass 3.3 GL 45 PU 10 pieces

**€35.30**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Graduated laboratory bottles with screw thread made of borosilicate glass 3.3 (acc. to ISO 4796 DIN 168 part 1).

## Additional Information

---

No.	2-3062
Manufacturer (Brand)	neoLab
EAN	4058072000486
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Volume	500 ml
Color bottle	transparent (clear)
Height	181 mm
Bottles properties	with high hydrolytic resistance with labeling field with pouring ring graduated
DM bottle	86 mm
DM Neck	45 mm
Thread	DIN 45
Neck size	narrow
Material bottle	Borosilicate glass (1st hydrolytic class)

