



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## DURAN® Laboratory bottle Duran GL 25, 25 ml with pouring ring and screw cap

**€90.60**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Original <sup>DURAN®</sup> quality laboratory bottles.

With DIN thread, clear glass. Complete with blue screw cap with integrated lip seal.

**Without pouring ring**

Uniform, strong wall. According to ISO 4796-1, autoclavable (open cap for this purpose). Caps are temperature resistant up to +140°C.

With batch identification "Retrace Code", batch certificate available online.

## Additional Information

No.	SD-0790
Manufacturer (Brand)	DURAN
old neoLab article no.	E-2055
EAN	4058072345273
Transport temperature	Room temperature
Color	transparent (clear)
Volume	25 ml
sterile	No
Color bottle	transparent (clear)
Height	74 mm
Type closure/ lid	Screw cap
Bottles properties	with high hydrolytic resistance with closure/lid with seal with batch certificate with labeling field graduated cap autoclavable bottle autoclavable
DM bottle	36 mm
DM Neck	25 mm
Color closure/ lid	blue
Thread	GL 25
Neck size	narrow
Material gasket / septum	Polypropylene (PP)
Material bottle	Borosilicate glass (1st hydrolytic class)
Material closure/ lid	Polypropylene (PP)
TBST MAX	525 °C
TBST MIN	-70 °C
Closure / lid properties	closed with seal

