



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## DURAN® Laboratory bottle GL 45, 250 ml with pouring ring and screw cap

**€12.04** As low as  
**€9.26**  
plus VAT &  
Shipping

### Product Images



## Description

---

Original DURAN® quality laboratory bottles.

With DIN thread, clear glass. Complete with blue screw cap with integrated lip seal and PP pouring ring.

Uniform, strong wall. According to ISO 4796-1, autoclavable (open cap for this purpose). Cap and pouring ring are temperature resistant up to +140°C.

With batch identification "Retrace Code", batch certificate available online.

## Additional Information

No.	SD-0800
Manufacturer (Brand)	DURAN
old neoLab article no.	E-2058
EAN	4058072345358
Transport temperature	Room temperature
Color	transparent (clear)
Volume	250 ml
sterile	No
Color bottle	transparent (clear)
Height	143 mm
Type closure/ lid	Screw cap
TBST MAX	525 °C
TBST MIN	-70 °C
Bottles properties	with high hydrolytic resistance with closure/lid with seal with batch certificate with labeling field graduated cap autoclavable bottle autoclavable
Thread	GL 45
Neck size	narrow
Material gasket / septum	Polypropylene (PP)
Closure / lid properties	closed with seal
DM bottle	70 mm
DM Neck	45 mm
Color closure/ lid	blue
Material bottle	Borosilicate glass (1st hydrolytic class)
Material closure/ lid	Polypropylene (PP)

