



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



DURAN® Laboratory bottle, clear, graduated, GL 45, with cap and pouring ring (PP), 2000 ml

€49.96 As low
as
€38.43
plus VAT &
Shipping

Product Images



Description

Original ^{DURAN®} quality laboratory bottles.

With DIN thread, clear glass. Complete with blue screw cap with integrated lip seal and PP pouring ring.

Uniform, strong wall. According to ISO 4796-1, autoclavable (open cap for this purpose). Cap and pouring ring are temperature resistant up to +140°C.

With batch identification "Retrace Code", batch certificate available online.

Additional Information

No.	SD-0811
Manufacturer (Brand)	DURAN
old neoLab article no.	E-2061
EAN	4058072345426
Transport temperature	Room temperature
Color	transparent (clear)
Volume	2000 ml
sterile	No
Color bottle	transparent (clear)
Height	265 mm
Type closure/ lid	Screw cap
TBST MAX	525 °C
TBST MIN	-70 °C
Bottles properties	with high hydrolytic resistance with closure/lid with seal with batch certificate with labeling field graduated cap autoclavable bottle autoclavable
Thread	GL 45
Neck size	narrow
Material gasket / septum	Polypropylene (PP)
Closure / lid properties	closed with seal
DM bottle	136 mm
DM Neck	45 mm
Color closure/ lid	blue
Material bottle	Borosilicate glass (1st hydrolytic class)
Material closure/ lid	Polypropylene (PP)

