



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e


Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**neolab® Laboratory bottles, borosilicate glass,  
GL 45, 1000 ml, cap + pouring ring, 10 pcs./pack**

**€72.70  
plus VAT &  
Shipping**

## Product Images



## Description

---

### Laboratory bottles

Simax®, borosilicate glass 3.3,  
DIN thread GL 45, autoclavable,  
Screw cap made of PPN,  
Pouring ring

## Additional Information

---

No.	E-1432
Manufacturer (Brand)	neoLab
EAN	4058072064723
Transport temperature	Room temperature
Color	transparent (clear)
Volume	1000 ml
sterile	No
Color bottle	transparent (clear)
Height	230 mm
Type closure/ lid	Screw cap
Bottles properties	with high hydrolytic resistance with cap/lid with labeling field with pouring ring graduated cap autoclavable bottle autoclavable
Thread	GL 45
Neck size	narrow
Closure / lid properties	closed
DM bottle	101 mm
DM Neck	45 mm
Color closure/ lid	blue
Material bottle	Borosilicate glass (1st hydrolytic class)
Material closure/ lid	Polypropylene (PP)

