



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e


Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® Laboratory bottles, borosilicate glass,
GL 45, 2000 ml, cap + pouring ring, 10 pcs./pack**

**€180.60
plus VAT &
Shipping**

Product Images



Description

Laboratory bottles

Simax®, borosilicate glass 3.3,
DIN thread GL 45, autoclavable,
Screw cap made of PPN,
Pouring ring

Additional Information

No.	E-1458
Manufacturer (Brand)	neoLab
EAN	4058072064921
Transport temperature	Room temperature
Color	transparent (clear)
Volume	2000 ml
sterile	No
Color bottle	transparent (clear)
Height	260 mm
Type closure/ lid	Screw cap
Bottles properties	with high hydrolytic resistance with cap/lid with labeling field with pouring ring graduated cap autoclavable bottle autoclavable
DM bottle	136 mm
DM Neck	45 mm
Color closure/ lid	blue
Thread	GL 45
Neck size	narrow
Material bottle	Borosilicate glass (1st hydrolytic class)
Material closure/ lid	Polypropylene (PP)
Closure / lid properties	closed

