



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Wheaton Media/Laboratory bottles with cap 219819 Cl.-Gl.500 ml Pack of 6 9-0195

€349.68
plus VAT &
Shipping

Product Images



Description

Laboratory media bottles made of Wheaton 400 borosilicate glass with low alkali content, which avoids changes in pH values and ensures the purity of the contents.

Graduation with maximum volume leaves enough space for mixing and shaking.

The white labeling field stands out from the bottle with high contrast and can be marked with many writing instruments.

With pouring lip.

The 1 l bottle has practical recessed grips.

Black screw caps with septum made of rubber/Teflon®.

Bottles and caps can be autoclaved.

Additional Information

No.	WD-0490
Manufacturer (Brand)	Wheaton
old neoLab article no.	9-0194
EAN	4058072426804
Transport temperature	Room temperature
Color	transparent (crystal clear)
Volume	500 ml
sterile	No
Color bottle	transparent (crystal clear)
Height	188 mm
Type closure/ lid	Screw cap
Bottles properties	with high hydrolytic resistance with closure/lid with recessed grip with labeling field with pouring ring graduated closure with septum closure autoclavable bottle autoclavable
Thread	33-430 (US)
Neck size	narrow
Material gasket / septum	Polytetrafluoroethylene (PTFE, Teflon®) butyl rubber (IIR)
Closure / lid properties	closed with mounted septum
DM bottle	88 mm
DM Neck	33 mm
Color closure/ lid	black
Material bottle	Borosilicate glass (Wheaton 400, USP standard type 1) Borosilicate glass (1st hydrolytic class)
Material closure/ lid	Phenolic resin (PF resin)

