



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Wheaton Media/Laboratory bottles with closure 1000 ml, pack of 6

€129.31 As low
as
€105.66
plus VAT &
Shipping



Description

Laboratory media bottles made of Wheaton 400 borosilicate glass with low alkali content, which avoids changes in pH values and ensures the purity of the contents.

Graduation with maximum volume leaves enough space for mixing and shaking.

The white labeling field stands out from the bottle with high contrast and can be marked with many writing instruments.

With pouring lip.

The 1 l bottle has practical recessed grips.

Black screw caps with septum made of rubber/Teflon®.

Bottles and caps can be autoclaved.

Additional Information

No.	WD-0491
Manufacturer (Brand)	Wheaton
old neoLab article no.	9-0195
EAN	4058072426811
Transport temperature	Room temperature
Color	transparent (crystal clear)
Volume	1000 ml
sterile	No
Color bottle	transparent (crystal clear)
Height	231 mm
Type closure/ lid	Screw cap
Bottles properties	with high hydrolytic resistance with closure/lid with recessed gripl with labeling field with pouring ring graduated closure with septum closure autoclavable bottle autoclavable
Thread	38-430 (US)
Neck size	narrow
Material gasket / septum	Polytetrafluoroethylene (PTFE, Teflon®) butyl rubber (IIR)
Closure / lid properties	closed with mounted septum
DM bottle	103 mm
DM Neck	38 mm
Color closure/ lid	black
Material bottle	Borosilicate glass (Wheaton 400, USP standard type 1) Borosilicate glass (1st hydrolytic class)
Material closure/ lid	Phenolic resin (PF resin)

