



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Rasotherm laboratory bottles, GL45, blue screw cap, (Boro 3.3), 500 ml

Special Price

**€4.00** was

**€5.18** As low as

**€3.58**

**plus VAT &  
Shipping**

## Product Images

---



## Description

---

The RASOTHERM® **laboratory thread bottles** GL 45 are made of particularly heat- & chemical-resistant borosilicate glass 3.3.

All bottles meet the highest quality standards and comply with ISO 4796-1.

Matching blue lids and PP pouring rings are included.

The bottles are graduated.

## Additional Information

No.	1-0197
Manufacturer (Brand)	Rasotherm
EAN	4058072081195
Color	transparent (clear)
Material	Borosilicate glass
Volume	500 ml
sterile	No
Color bottle	transparent (clear)
Height	182 mm
Type closure/ lid	Screw cap
Bottles properties	with closure/lid with pouring ring graduated closure autoclavable bottle autoclavable
Material	Borosilicate glass
Thread	GL 45
Neck size	narrow
Closure / lid properties	closed
DM bottle	86 mm
DM Neck	45 mm
Spare part	1-0202
Color closure/ lid	blue
Material bottle	Borosilicate glass
Material closure/ lid	Polypropylene (PP)

