



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Rasotherm laboratory bottles, (Boro 3.3), GL45, blue screw cap, 1000 ml

Special Price

€6.50 was

€8.27 As low as

€5.94

**plus VAT &
Shipping**

Product Images



Description

The RASOTHERM® **laboratory thread bottles** GL 45 are made of particularly heat- & chemical-resistant borosilicate glass 3.3.

All bottles meet the highest quality standards and comply with ISO 4796-1.

Matching blue lids and PP pouring rings are included.

The bottles are graduated.

Additional Information

No.	1-0198
Manufacturer (Brand)	Rasotherm
EAN	4058072080761
autoclavable at 121°C	No
Color	transparent (clear)
Material	Borosilicate glass
Volume	1000 ml
sterile	No
Color bottle	transparent (clear)
Height	230 mm
Type closure/ lid	Screw cap
Bottles properties	with closure/lid with pouring ring graduated closure autoclavable bottle autoclavable
Material	Borosilicate glass
Thread	GL 45
Neck size	narrow
Closure / lid properties	closed
DM bottle	101 mm
DM Neck	45 mm
Color closure/ lid	blue
Spare part	1-0202
Material bottle	Borosilicate glass
Material closure/ lid	Polypropylene (PP)

