



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Rasotherm laboratory bottles, (Boro 3.3), GL45, blue screw cap, 5000 ml

Special Price

€65.30 was

€83.70

**plus VAT &
Shipping**



Description

The RASOTHERM® **laboratory thread bottles** GL 45 are made of particularly heat- & chemical-resistant borosilicate glass 3.3.

All bottles meet the highest quality standards and comply with ISO 4796-1.

Matching blue lids and PP pouring rings are included.

The bottles are graduated.

Additional Information

| | |
|--------------------------|---|
| No. | 1-0200 |
| Manufacturer (Brand) | Rasotherm |
| EAN | 4058072080785 |
| Color | transparent (clear) |
| Material | Borosilicate glass |
| Volume | 5000 ml |
| sterile | No |
| Color bottle | transparent (clear) |
| Height | 330 mm |
| Type closure/ lid | Screw cap |
| Bottles properties | with closure/lid with pouring ring graduated closure autoclavable bottle autoclavable |
| Material | Borosilicate glass |
| Thread | GL 45 |
| Neck size | narrow |
| Closure / lid properties | closed |
| DM bottle | 186 mm |
| DM Neck | 45 mm |
| Spare part | 1-0202 |
| Color closure/ lid | blue |
| Material bottle | Borosilicate glass |
| Material closure/ lid | Polypropylene (PP) |

