



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Socorex PP container, 250 ml

€36.61
plus VAT &
Shipping

Product Images

Description

PP reservoirs for Socorex Dosys™ self-filling syringes

Large screw-on PP reservoirs with air supply are available in three different sizes.

Can only be used with Rotlauf aspiration valve (to be ordered separately under part no. 7-7711)

Autoclavable.

Additional Information

No.	SR-0307
Manufacturer (Brand)	Socorex
old neolab article no.	7-7710
EAN	4058072368241
Transport temperature	Room temperature
Color	transparent (milky)
Volume	250 ml
Bottles properties	autoclavable
Material	Polypropylene (PP)
Floor outside	rounded transition to the bottle bottom (HS bottom)
Floor inside	rounded transition to the bottle bottom (HS bottom)

