



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**Lenz oxygen cylinder acc. to Winkler, nominal  
volume 100 - 150 ml, ground joint NS 14/23**

**€27.45  
plus VAT &  
Shipping**

## Product Images



## Description

---

### Oxygen bottles acc. to Winkler

For the determination of the oxygen dissolved in water; content engraving on the bottle (accuracy:  $\pm 0.01$  ml); with labeling field; with solid, obliquely cut NS glass stopper (can be secured with spring clip); bottle and corresponding stopper adjusted together and labeled with the same number

## Additional Information

---

No.	LZ-1556
Manufacturer (Brand)	Lenz
old neoLab article no.	B-3559
EAN	4058072312947
Transport temperature	Room temperature
Color	transparent (clear)
Volume	150 ml
Color bottle	transparent (clear)
Type closure/ lid	NS plug
Bottles properties	with cap/lid cap autoclavable bottle autoclavable
Neck size	narrow
Color closure/ lid	transparent (clear)
Material bottle	Lime soda glass
Material closure/ lid	Glass

