



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**DURAN® Square screw cap bottle, wide neck,
clear, 250 ml, without pouring ring and dust
cap, soda lime glass**

**€55.10
plus VAT &
Shipping**

Product Images



Description

Square screw top bottle, soda-lime glass, clear, wide neck, low form, with DIN 45 S thread, 250 ml square bottle made of soda-lime clear glass (type 3) with wide neck. Special features and advantages: Ergonomic handling, high stability and good stackability due to angular shape Easy access to the bottle contents due to the wide diameter of the bottle neck Application: storage of samples and reagents, sampling, transport and shipping

Additional Information

No.	SD-1088
Reg. trademark	DURAN ®
Manufacturer (Brand)	DURAN
VGKL number	122382036
EAN	4058072188207
Manufacturer EAN	4032051024854
Gross weight	0.303
Storage temperature	Room temperature
Transport temperature	Room temperature
Net weight	0.28000 kg
Color	Crystal clear
Material	Lime soda glass
Volume	250 ml
Country of origin	Czech Republic
sterile	No
Packaging volume	0.00913300 cc
Customs tariff number	70179000
Neck size	far

