



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



DURAN® Square screw top vials amber glass 100 ml narrow neck, 10 pcs./pack

€45.90
plus VAT &
Shipping

Product Images



Description

Amber glass narrow-neck bottles with DIN thread.

Made of soda-lime glass of hydrolytic class 3.

Thermal shock resistance 30 K, a maximum service temperature of +100°C is recommended.

Please order closures, pouring rings, etc. separately.

From the DURAN® Group.

Additional Information

No.	SD-1082
Manufacturer (Brand)	DURAN
old neolab article no.	2-1375
EAN	4058072347840
Transport temperature	Room temperature
Color	brown
Material	Soda-lime glass (normal glass) 2
Volume	100 ml
Wide	49 mm
Color bottle	brown
Height	119 mm
Depth	49 mm
Bottles properties	with batch certificate bottle autoclavable
Material	Soda-lime glass (normal glass) 2
DM Neck	32 mm
Thread	DIN 32
Neck size	narrow
Material bottle	Soda-lime glass (normal glass) 2

