



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



DURAN® Square screw top vials amber glass 250 ml narrow neck, 10 pcs./pack

€10.36 As low as
€7.97
plus VAT &
Shipping

Product Images



Description

Amber glass narrow-neck bottles with DIN thread.

Made of soda-lime glass of hydrolytic class 3.

Thermal shock resistance 30 K, a maximum service temperature of +100°C is recommended.

Please order closures, pouring rings, etc. separately.

From the DURAN® Group.

Additional Information

| | |
|------------------------|--|
| No. | SD-1083 |
| Manufacturer (Brand) | DURAN |
| old neolab article no. | 2-1376 |
| EAN | 4058072347857 |
| Transport temperature | Room temperature |
| Color | brown |
| Material | Soda-lime glass (normal glass) 2 |
| Volume | 250 ml |
| Wide | 64 mm |
| Color bottle | brown |
| Height | 155 mm |
| Depth | 64 mm |
| Bottles properties | with batch certificate bottle autoclavable |
| Material | Soda-lime glass (normal glass) 2 |
| Thread | DIN 32 |
| Neck size | narrow |
| DM Neck | 32 mm |
| Material bottle | Soda-lime glass (normal glass) 2 |

