



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



DURAN® Square screw top bottles amber glass 500 ml narrow neck, 10 pcs./pack

€11.49 As low as
€8.84
plus VAT &
Shipping

Product Images



Description

Amber glass narrow-neck bottles with DIN thread.

Made of soda-lime glass of hydrolytic class 3.

Thermal shock resistance 30 K, a maximum service temperature of +100°C is recommended.

Please order closures, pouring rings, etc. separately.

From the DURAN® Group.

Additional Information

No.	SD-1084
Manufacturer (Brand)	DURAN
old neoLab article no.	2-1377
EAN	4058072347864
Transport temperature	Room temperature
Color	brown
Material	Soda-lime glass (normal glass) 2
Volume	500 ml
Wide	77 mm
Color bottle	brown
Height	186 mm
Depth	77 mm
Bottles properties	with batch certificate bottle autoclavable
Material	Soda-lime glass (normal glass) 2
Thread	DIN 32
Neck size	narrow
DM Neck	32 mm
Material bottle	Soda-lime glass (normal glass) 2

