



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## DURAN® Square screw top bottles amber glass 1000 ml narrow neck, 10 pcs./pack

**€17.45** As low as  
**€13.42**  
plus VAT &  
Shipping

### Product Images



## Description

---

Amber glass narrow-neck bottles with DIN thread.

Made of soda-lime glass of hydrolytic class 3.

Thermal shock resistance 30 K, a maximum service temperature of +100°C is recommended.

Please order closures, pouring rings, etc. separately.

From the DURAN® Group.

## Additional Information

---

No.	SD-1085
Manufacturer (Brand)	DURAN
old neolab article no.	2-1378
EAN	4058072347871
Transport temperature	Room temperature
Color	brown
Material	Soda-lime glass (normal glass) 2
Volume	1000 ml
Wide	97 mm
Color bottle	brown
Height	233 mm
Depth	97 mm
Bottles properties	with batch certificate bottle autoclavable
Material	Soda-lime glass (normal glass) 2
Thread	DIN 45
Neck size	narrow
DM Neck	45 mm
Material bottle	Soda-lime glass (normal glass) 2

