



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## DURAN® Square screw top vials clear glass 100 ml narrow neck, 10 pcs./pack

**€45.90**  
plus VAT &  
Shipping

### Product Images



## Description

---

### Narrow neck bottles with DIN thread, original DURAN®

Made of soda-lime glass, clear, hydrolytic class 3.

Thermal shock resistance 30 K, a maximum service temperature of +100°C is recommended.

Please order closures, pouring rings etc. separately.

## Additional Information

---

No.	SD-1078
Manufacturer (Brand)	DURAN
old neolab article no.	2-1371
EAN	4058072347802
Transport temperature	Room temperature
Color	transparent (clear)
Material	Soda-lime glass (normal glass) 2
Volume	100 ml
Wide	49 mm
Color bottle	transparent (clear)
Height	119 mm
Depth	49 mm
Bottles properties	with batch certificate bottle autoclavable
Material	Soda-lime glass (normal glass) 2
DM Neck	32 mm
Thread	DIN 32
Neck size	narrow
Material bottle	Soda-lime glass (normal glass) 2

