



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



DURAN® Square screw top vials clear glass 250 ml narrow neck, 10 pcs./pack

€79.70
plus VAT &
Shipping

Product Images



Description

Narrow neck bottles with DIN thread, original DURAN®

Made of soda-lime glass, clear, hydrolytic class 3.

Thermal shock resistance 30 K, a maximum service temperature of +100°C is recommended.

Please order closures, pouring rings etc. separately.

Additional Information

No.	SD-1079
Manufacturer (Brand)	DURAN
old neolab article no.	2-1372
EAN	4058072347819
Transport temperature	Room temperature
Color	transparent (clear)
Material	Soda-lime glass (normal glass) 2
Volume	250 ml
Wide	64 mm
Color bottle	transparent (clear)
Height	155 mm
Depth	64 mm
Bottles properties	with batch certificate bottle autoclavable
Material	Soda-lime glass (normal glass) 2
Thread	DIN 32
Neck size	narrow
DM Neck	32 mm
Material bottle	Soda-lime glass (normal glass) 2

