



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Wheaton square bottles 125 ml PP screw cap, tef.insert, VE24 Nr9-0072

€60.24
plus VAT &
Shipping

Product Images



Description

Bottles made of soda-lime glass type III, chemically particularly resistant.

Includes screw cap with PE foam insert and Teflon® coating.

Space-saving due to square shape.

Additional Information

No.	WD-1805
Manufacturer (Brand)	Wheaton
old neoLab article no.	9-0071
EAN	4058072427894
Transport temperature	Room temperature
Color	transparent (clear)
Volume	125 ml
Wide	45 mm
Color bottle	transparent (clear)
Height	111 mm
Depth	45 mm
Bottles properties	with closure/lid closure with septum bottle autoclavable
Thread	33-400 (US)
Neck size	narrow
Material gasket / septum	Polytetrafluoroethylene (PTFE, Teflon®) High-density polyethylene (HDPE)
Closure / lid properties	closed
DM Neck	33 mm
Color closure/ lid	white
Material bottle	Soda-lime glass (Wheaton 800, USP standard type 3)
Material closure/ lid	Polypropylene (PP)

