



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



## Wheaton square bottles 500 ml PP screw cap, tef.insert, VE24 Nr9-0074

**€104.88**  
**plus VAT &**  
**Shipping**

### Product Images



## Description

---

Bottles made of soda-lime glass type III, chemically particularly resistant.

Includes screw cap with PE foam insert and Teflon® coating.

Space-saving due to square shape.

## Additional Information

---

No.	WD-1807
Manufacturer (Brand)	Wheaton
old neoLab article no.	9-0073
EAN	4058072427917
Transport temperature	Room temperature
Color	transparent (clear)
Volume	500 ml
Wide	68 mm
Color bottle	transparent (clear)
Height	167 mm
Depth	68 mm
Bottles properties	with closure/lid closure with septum bottle autoclavable
DM Neck	43 mm
Color closure/ lid	white
Thread	48-400 (US)
Neck size	narrow
Material gasket / septum	Polytetrafluoroethylene (PTFE, Teflon®) High-density polyethylene (HDPE)
Material bottle	Soda-lime glass (Wheaton 800, USP standard type 3)
Material closure/ lid	Polypropylene (PP)
Closure / lid properties	closed

