




neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



Falcon® storage bottles for cell culture media, 500 ml, 24 pcs./pack

€151.27
plus VAT &
Shipping

Product Images



Description

Disposable polystyrene bottles, ideal for storing media, buffer solutions and other aqueous solutions. Indented grip surfaces facilitate handling. The airtight lid minimizes the risk of contamination. The bottles can be used excellently in conjunction with the Bottle Top Filters. They are sterile, certified and pyrogen-free.

Additional Information

No.	XC-0691
Manufacturer (Brand)	Falcon
old neoLab article no.	7-4182
EAN	4058072444167
Transport temperature	Room temperature
Color	transparent (crystal clear)
Material	Polystyrene (PS)
Volume	500 ml
sterile	Yes
Color bottle	transparent (crystal clear)
Type closure/ lid	Screw cap
Bottles properties	with closure/lid with recessed grip graduated closure autoclavable
Material	Polystyrene (PS)
Neck size	narrow
Closure / lid properties	closed
DM Neck	45 mm
Color closure/ lid	orange
Material bottle	Polystyrene (PS)
Material closure/ lid	Polypropylene (PP)

