



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**DURAN® Wide-neck stand-up bottle, NS 60/46,
with stopper, clear, soda-lime glass, 1000 ml**

**€366.30
plus VAT &
Shipping**

Product Images



Description

Stand-up bottle, wide neck, soda-lime glass, clear, with standard ground joint, with glass stopper, 1000 ml A maximum usage temperature of +100 °C is recommended. Thermal shock resistance 30 K. Hydrolytic class 3. Note on using DURAN® glass stoppers with soda-lime glass bottles: If the bottle and the stopper have a temperature difference greater than 30 °C, the glass stoppers can become stuck! Special features and advantages: Easy access to the bottle contents due to the wide diameter of the bottle neck Consists only of glass components Tight closure due to standard ground joint Application: Storage of air-sensitive reagents and chemicals, storage of hygroscopic liquids, granulates and powders

Additional Information

No.	SD-0964
Reg. trademark	DURAN ®
Manufacturer (Brand)	DURAN
VGKL number	122318554
EAN	4058072174873
Manufacturer EAN	4032051022751
Gross weight	0.83
Storage temperature	Room temperature
Transport temperature	Room temperature
Net weight	0.77500 kg
Color	Crystal clear
Material	Lime soda glass
Volume	1000 ml
Country of origin	Germany
sterile	No
Region of origin	Baden-Württemberg
Packaging volume	0.02948100 cc
Customs tariff number	70172000
Neck size	far
Standard ground	60/46

