



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



DURAN® Laboratory bottles amber glass GLS 80, 1000 ml without pouring ring and screw cap

€89.58 As low
as
€68.91
plus VAT &
Shipping

Product Images



Description

DURAN®^{laboratory} bottle with wide thread GLS 80 (= 80 mm OD) made of amber glass, based on ISO 4796.

A special thread allows opening with less than one turn. The wide neck allows convenient filling and removal of powders and viscous substances.

The brown coloration provides light protection from approx. 300-500 nm. Unchanged DURAN® properties on the inside, as color application is only on the outside. Long shelf life and good chemical resistance of the brown coloring.

Large base area for safe standing and stackability.

White text field, easy-to-read graduation Batch identification "Retrace Code" Batch certificate available online. Inert and chemically resistant.

Accessories: Screw cap (Order No. 2-3059)

Pouring ring (Order No. 2-3049)

Additional Information

No.	SD-0130
Manufacturer (Brand)	DURAN
old neoLab article no.	E-2041
EAN	4058072339937
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	brown
Volume	1000 ml
Color bottle	brown
Height	248 mm
Bottles properties	with high hydrolytic resistance with batch certificate with labeling field graduated
Thread	GLS 80
Neck size	far
DM bottle	136 mm
DM Neck	80 mm
Material bottle	Borosilicate glass (1st hydrolytic class)

