



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



DURAN® GLS 80 wide-neck bottle, clear, 500 ml, without cap and ring

€18.47 As low as
€14.21
plus VAT &
Shipping

Product Images



Description

A specially developed thread allows opening with less than one turn. The 80 mm wide outer diameter of the bottle neck allows convenient, safe filling and removal of powders and viscous substances. Even large-volume funnels can be fitted without difficulty.

Sturdy design gives high mechanical load capacity.

Large base area for safe standing and stackability.

White lettering field, easy-to-read graduation, batch identification "Retrace Code". Inertness and chemical resistance due to DURAN®

Accessories: Screw cap (Order No. 2-3059)

Pouring ring (Order No. 2-3049)

Additional Information

No.	SD-0123
Manufacturer (Brand)	DURAN
old neoLab article no.	2-3053
EAN	4058072339869
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	transparent (clear)
Volume	500 ml
Color bottle	transparent (clear)
Height	148 mm
TBST MAX	525 °C
TBST MIN	-70 °C
Bottles properties	with high hydrolytic resistance with batch certificate with labeling field graduated
Thread	GLS 80
Neck size	far
DM bottle	101 mm
DM Neck	80 mm
Material bottle	Borosilicate glass (1st hydrolytic class)

