



neolab Migge GmbH  
Rischerstr. 7-9  
69123 Heidelberg  
Deutschland  
+49 (0)6221 /  
8442-44  
<https://www.neolab.de>  
e

Umsatzsteuer-  
Identifikationsnummer  
:  
DE 143 450 657



**neolab® Wide-neck screw cap bottle 100 ml,  
amber glass, DIN 45, without screw cap**

**€1.80** As low as  
**€1.47**  
**plus VAT &  
Shipping**

## Product Images

---



## Description

---

Wide neck packaging bottles made of brown soda glass, round, according to DIN.  
Accessories: Screw thread caps according to DIN (Order No. 2-1495 ff)

## Additional Information

---

No.	DT-0083
Manufacturer (Brand)	neoLab
old neoLab article no.	2-1481
EAN	4058072228590
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	brown
Material	Soda-lime glass (normal glass) 2
Volume	100 ml
Color bottle	brown
Height	85.5 mm
Material	Soda-lime glass (normal glass) 2
Thread	DIN 45
Neck size	far
DM bottle	50 mm
DM Neck	45 mm
Material bottle	Soda-lime glass (normal glass) 2

