



neolab Migge GmbH
Rischerstr. 7-9
69123 Heidelberg
Deutschland
+49 (0)6221 /
8442-44
<https://www.neolab.de>
e

Umsatzsteuer-
Identifikationsnummer
:
DE 143 450 657



**neolab® Wide-neck screw cap bottle 250 ml,
amber glass, DIN 55, without screw cap**

€2.66 As low as

€2.28

**plus VAT &
Shipping**

Product Images



Description

Wide neck packaging bottles made of brown soda glass, round, according to DIN.
Accessories: Screw thread caps according to DIN (Order No. 2-1495 ff)

Additional Information

No.	DT-0084
Manufacturer (Brand)	neoLab
old neoLab article no.	2-1482
EAN	4058072228606
Transport temperature	Room temperature
autoclavable at 121°C	Yes
Color	brown
Material	Soda-lime glass (normal glass) 2
Volume	250 ml
Color bottle	brown
Height	113 mm
Material	Soda-lime glass (normal glass) 2
DM bottle	63 mm
DM Neck	55 mm
Thread	DIN 55
Neck size	far
Material bottle	Soda-lime glass (normal glass) 2

